Saypha® VOLUME

EN

Instructions for use

Composition:
Highly cross-linked sodium hyaluronate 23 mg/mL. Phosphate buffer pH 6.8–7.4 q.s., sodium chloride. One syringe contains 1.0 mL of Saypha® VOLUME. Sterilized using steam.

Description:
Saypha® VOLUME is a sterile, biodegradable, viscoelastic, clear, transparent, isotonic and homogenized, injectable gel implant. Saypha® VOLUME consists of cross-linked hyaluronic acid (HA), obtained from *Streptococcus equi* bacteria, formulated to a concentration of 23 mg/mL in a physiologic buffer. Each box contains one syringe of Saypha® VOLUME, two 27G ½" disposable sterile needles reserved for injection of Saypha® VOLUME and a product leaflet. A set of two labels showing the batch number is situated at the bottom of the box. One of these labels should be attached to the patient's file and the other should be given to the patient to ensure traceability.

Indications:
Saypha® VOLUME is a viscoelastic solution to correct deeper wrinkles in nasolabial folds. It is indicated to be injected into the deep dermis.

Exclusion criteria:
Do not inject Saypha® VOLUME into the periorbital area (crow's feet, eye circle, eye lids) or glabellar region. Saypha® VOLUME must not be used in:
- patients who tend to develop hypertrophic scarring. | - patients with a history of autoimmune disease or who are receiving immune therapy. | - patients who are known to be hypersensitive to hyaluronic acid. | - pregnant or breastfeeding women. | - patients under 18 years of age.
Patients on anticoagulants or patients receiving platelet aggregation inhibitors (e.g. ASS) should not be treated with Saypha® VOLUME without consulting their doctors.
Saypha® VOLUME must not be used in areas presenting cutaneous, inflammatory and/or infectious processes (e.g. acne, herpes ...). Saypha® VOLUME must not be used in association with laser therapy, chemical peeling or dermabrasion.

Precautions to use:
Saypha® VOLUME is only indicated for intradermal injections. Sensitive skin may be pretreated using a local anaesthetic patch or cream. Please note that any anaesthesia may cause redness or local hypersensitivity. There are no available clinical data (efficiency, tolerance) about injecting Saypha® VOLUME into an area which has already been treated with another filling product. Patients should be advised not to apply any make-up for 12 hours after the injection and to avoid prolonged exposure to sunlight and UV or using saunas or Turkish baths for one week after the injection. If the 27G ½" needle is blocked, do not increase the pressure on the plunger rod but stop the injection and replace the needle. Do not inject into blood vessels. Do not use Saypha® VOLUME in bones, tendons, ligaments or muscles. Do not inject Saypha® VOLUME into naevi. Do not overcorrect. Discard the syringe and remaining product after use. Do not use if packaging has been damaged.

Interactions with other agents:
There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore Saypha® VOLUME should never be placed in contact with these substances or with medical-surgical instruments that have been in contact with these substances.

Undesired side effects:
Physicians must inform the patient that there are potential side effects and/or incompatibilities associated with the implantation of this device, which may occur immediately or may be delayed. These include (non-exhaustive list):
- Inflammatory reactions (e.g. redness, oedema, erythema ...) which may be associated with itching and pain on pressure after the injection. These reactions may last for a week. | - Indurations or nodules at the injection site especially in cases of too superficial placement. In a clinical study of similar substances, local mobility of the injected gel and skin limpness were observed; this might be due to the use of too much volume and/or sub-optimal injection technique (e.g. large bolus injections). In some cases implant aspiration was used to resolve the event. | - Discolouration of the injection site. | - Poor effect or weak filling effect if Saypha® VOLUME is injected incorrectly.
Cases of skin necrosis, bacterial infection, formation of granulomas, cystic and acneiform lesions and hypersensitivity have all been reported in the literature following hyaluronic acid injections. Cases of visual loss as the most severe uncommon undesired events due to wrong injection technique have also been reported. It is therefore important to take these possible complications into account. Patients must inform their physician as soon as possible about any inflammatory reactions which persist for more than one week or any other secondary effect which develops. The physician should treat these as appropriate. Any other undesirable side effects associated with the injection of Saypha® VOLUME must be reported to the distributor and/or to the manufacturer.

Methods of use:
Saypha® VOLUME must be injected into non-inflamed, disinfected, healthy skin. The technique used is essential for the success of the treatment (e.g. microdroplet technique). Therefore this device must only be used by doctors who have received specific training on the injection technique to restore volume to the face. The area to be treated should be disinfected thoroughly prior to the injection. Use the 27G ½" needle which is provided with the syringe and inject slowly by applying the appropriate injection technique. The amount injected will depend on the wrinkle or depression which is to be corrected. A maximum volume of 2.0mL and 1.0mL per nasolabial fold is recommended for initial and touch-up treatments, respectively. After the injection, doctors may apply a light massage in order to distribute the product uniformly. However, patients should be advised not to massage by themselves for one week following the injection in order to avoid a possible risk of product mobility.

Warnings:
Verify the integrity of the syringe and the expiry date before use. Do not use any other needle or syringe than provided by the manufacturer. Do not re-use; quality and sterility can only be guaranteed for an originally closed syringe. The re-use of the product creates a potential infection risk for patients or users. Do not re-sterilize. After use the needle must be disposed of in a special container.

Storage:
Store Saypha® VOLUME at 2–25 °C / 36–77 °F in a dry place in the original box. Protect from light, heat and frost and handle with care.

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
Saypha® VOLUME

Instructions for the correct handling of the device:
Hold the Luer-Lock adapter as shown in ❶. To remove the tip cap, twist ❷ and pull carefully ❸. Following the instruction above will prevent the formation of air bubbles. Hold the syringe as shown in ❹. Open the enclosed needle case and insert the needle firmly ❺ (do not use any other needle). Tighten the cap by twisting it clockwise ❻.

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